



# TRAPPING FLOWPATH - ST CONFIGURATION

## LOW PRESSURE SELECTORS


### Technical Note 713




**FIGURE 1: Valve body**



**FIGURE 2: As installed on Electric Actuators**



**FIGURE 2: As installed on Air Actuators**



CST (1/16") and ST (1/8") selectors are used for multi-column, multi-sample, or multi-trap operations, and are available for use with 4 to 16 loops, or positions.

Devices to be selected (loops, columns, spargers, etc. are installed in each pair of ports. Refer to the drawing above to determine which ports are paired.

**Watch an animation of the ST selector flowpath at [vici.com](http://vici.com):**



## SAMPLE TRAPPING APPLICATIONS FOR 4 TO 16 STREAMS

ST selectors are used for multi-column, multi-sample, or multi-trap operations. The ST configuration is available in both MW and UW type designs.

A typical application, shown here, is the collection of fractions at timed intervals for analysis at a later time. Valves can be ordered with matched loops already installed.

In this example, the 6 port valve shown is used to select between **1** collection/trapping and **2** analysis/desorption. Both valves can be supplied with pneumatic or electric actuators to automate these functions.

Watch an animation of the ST selector flowpath at [vici.com](http://vici.com):



### ADDITIONAL TECHNICAL NOTES

Technical notes related to cleaning and rotor replacement and to proper valve/actuator alignment can be found in the Support section of [vici.com](http://vici.com).

### TECH TIP

Because the most common cause of valve failure is stray particulates entering the valve, we strongly recommend the use of in-line filters at sample entry points.

Our ZUFR filters feature inexpensive and easily replaceable low pressure drop filter screens (2 or 10 micron). The filters are available in 1/16", 1/8", and 1/4" standard, reducing, and bulkhead versions.

## APPLICATION - ST FLOWPATH

