

CHEMINERT VALVE	X		Y	
	mm	Inch	mm	Inch
C5-* 4 or 6 port	ø31.7	ø1.25"	33	1.3"
C5H-* 8 or 10 port	ø34.9	ø1.375"	35.1	1.38"
C75-* 4 or 6 port	ø31.7	ø1.25"	38.1	1.5"
C75-* 8 or 10 port	ø34.9	ø1.375"	40.2	1.58"

MOTOR DRIVER CABLE  
PART NO. I-22640-CE  
3 FT LONG

POWER CABLE 24VDC

MOUNTING BRACKET  
PART NO. I-22766

MANUAL CONTROLLER  
PART NO. I-22679



DIMENSIONS:  
58 x 92 x 27 mm  
2.27 x 3.62 x 1.06 inch

**VICI** Valco Instruments Co. Inc.

MULTIPOSITION ACTUATOR  
CONTROL MODULE

CONTROL MODULE  
BOTTOM SIDE

REMOTE CONTROL CABLE  
PART NO. I-22696  
6 FT LONG



MANUAL CONTROLLER CABLE  
PART NO. I-22696-01  
3 FT LONG

CONTROL MODULE  
PART NO. FOR EMH (High speed): EMHCA-CE  
PART NO. FOR EMT (High torque): EMTCA-CE

DIMENSIONS:  
94 x 114 x 38.1 mm  
3.7 x 4.5 x 1.5 inch

Dimension in mm [Inch]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen: 1) 18.09.98 FP #90, chg. bracket part no. 2) 12.11.09 FP #602, chg. part no. valve 3) 25.03.13 FP #716, chg. Mem. controller design	Model: MTEVAL	Dimensions: ORIGINAL, YOU REQUEST
PART NO.: 18.11.1997 F.P. 25.03.2013 FCo 506011.DWG	Date: 18.11.1997 F.P. 25.03.2013 FCo 506011.DWG	Customer: KUNDE
<b>Micro electric actuator Multiposition valve</b> with close-mount hdw. General dimensions		Zeichnungs-Nr. / DRAWING NO. <b>1-GE-5060.11 d</b>

VICI® VICI AG International, CH-6214 Scherikon