



SHOWN ZP1



SHOWN ZP1H



\* different materials

Part#	Description	L		X	
		mm	Inch	mm	Inch
ZP.5	Zero volume plug, ZDV: 1/32"	~12.5	~0.49"	4.76	3/16" Hex
ZP.5H	Zero volume plug, ZDV: 1/32", HI-PSI	~12.5	~0.49"	4.76	3/16" Hex
ZP1*	Zero volume plug, ZDV: 1/16"	~19.1	~0.75"	6.35	1/4" Hex
ZP1H*	Zero volume plug, ZDV: 1/16", HI-PSI	~17.5	~0.69"	6.35	1/4" Hex
LZP1*	Zero volume plug, ZDV: 1/16", long	~28.8	~1.13"	6.35	1/4" Hex
LZP1H*	Zero vol. plug, ZDV: 1/16", long, HI-PSI	~25.7	~1.01"	6.35	1/4" Hex
ZP2*	Zero volume plug, ZDV: 1/8"	~25.4	~1.00"	9.53	3/8" Hex
ZP2H*	Zero volume plug, ZDV: 1/8", HI-PSI	~24.4	~0.96"	9.53	3/8" Hex
LZP2*	Zero volume plug, ZDV: 1/8", long	~35.6	~1.21"	9.53	3/8" Hex
ZP4	Zero volume plug, ZDV: 1/4"	~30.5	~1.20"	12.7	1/2" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material  
MATERIAL


Dimensionen  
DIMENSIONS

Aenderungen: a) 07.05.09 FP #577, added Part list  
MODIFICATIONS: b) 14.05.09 FP #579, added Part list

Projektion  
PROJECTION

Allgemein Toleranzen  
GENERAL TOLERANCES

PART NO.:



SN EN 22768-1/ ISO 2768

# Plugs, stainless steel and high pressure

General dimensions

Masstab  
SCALE

3:1

Gezeichnet  
DRAWN

21.12.2005

F.P.

Geprüft / Freigabe  
CHECKED/ RELEASE

14.05.2009

FCo

CAD Nr.  
CAD NO.

506065.idw

Unterverz.:  
SUB DIREC.:

KUNDE

Zeichnungs Nr. / DRAWING NO.

4 - GE - 5060.65 b



VICI AG International, CH-6214 Schenkon