


SHOWN CEF1.5



Part#	Description	Bore		L		X1	X2	X3
		mm	Inch	mm	Inch			
CEF1.5*	Internal column end fitting 1/16" to 1/32"	Ø0.25	Ø.010"	-32.0	-1.26"	1/4" Hex	1/4" Hex	3/16" Hex
CEF1*	Internal column end fitting 1/16" to 1/16"	Ø0.25	Ø.010"	-37.5	-1.48"	1/4" Hex	1/4" Hex	1/4" Hex
CEF2.51.0	Internal column end fitting 1/8" to 1/32"	Ø0.25	Ø.010"	-38.0	-1.50"	7/16" Hex	3/8" Hex	3/16" Hex
CEF211.0*	Internal column end fitting 1/8" to 1/16"	Ø0.25	Ø.010"	-43.5	-1.71"	7/16" Hex	3/8" Hex	1/4" Hex
CEF212.0*	Internal column end fitting 1/8" to 1/16"	Ø0.25	Ø.010"	-43.5	-1.71"	7/16" Hex	3/8" Hex	1/4" Hex
CEF412.1	Internal column end fitting 1/4" to 1/16"	Ø0.40	Ø.016"	-49.0	-1.93"	9/16" Hex	1/2" Hex	1/4" Hex
CEF413.*	Internal column end fitting 1/4" to 1/16"	Ø0.40	Ø.016"	-49.0	-1.93"	9/16" Hex	1/2" Hex	1/4" Hex
CEF414.*	Internal column end fitting 1/4" to 1/16"	Ø0.40	Ø.016"	-49.0	-1.93"	9/16" Hex	1/2" Hex	1/4" Hex
CEF616.0*	Internal column end fitting 3/8" to 1/16"	Ø0.51	Ø.020"	-53.5	-2.10"	3/4" Hex	11/16" Hex	1/4" Hex
CEF617.0*	Internal column end fitting 3/8" to 1/16"	Ø0.51	Ø.020"	-53.5	-2.10"	3/4" Hex	11/16" Hex	1/4" Hex
CEF819.0*	Internal column end fitting 1/2" to 1/16"	Ø0.76	Ø.030"	-56.0	-2.20"	7/8" Hex	13/16" Hex	1/4" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Part No. CEF*	Material PRODUCTION	Dimensional TOLERANCES
Modifications	Position	General Tolerances
Material	Scale	EN ISO 2768-V

Internal column end fitting		Material	11.02.2008	F.P.
Scale	6:1	Revision	11.02.2008	FCO
CAD No.	506319.idw	Author		
Customer	KUNDE	Drawn		