



Dimension in mm [Inch]



SHOWN EUTA-CSC6MWE

Type	Connections	No. of Pos.	X		Y		L	
			mm	Inch	mm	Inch	mm	Inch
CSD*MW	1/16"	4, 6, 8, 10, 12, 16	Ø38.1	Ø1.5"	~34.5	~1.36"	41.3	1.625"
SD*MW	1/8"	4, 6, 8, 10, 12	Ø50.8	Ø2.0"	~40.5	~1.60"	41.9	1.65"
SD*MW	1/8"	16	Ø57.15	Ø2.25"	~51.0	~2.00"	41.9	1.65"
CSC*MW	1/16"	4, 6, 8, 10, 12, 16	Ø41.3	Ø1.625"	~36.0	~1.42"	51.4	2.025"
SC*MW	1/8"	4, 6, 8, 10	Ø50.8	Ø2.0"	~38.5	~1.52"	52.1	2.05"
SC*MW	1/8"	12, 16	Ø57.15	Ø2.25"	~43.0	~1.69"	52.1	2.05"
CSF*MW	1/16"	4, 6, 8, 10, 12, 16	Ø41.3	Ø1.625"	~35.0	~1.38"	51.4	2.025"
SF*MW	1/8"	4, 6, 8, 10	Ø50.8	Ø2.0"	~39.5	~1.56"	52.1	2.05"
SF*MW	1/8"	12, 16	Ø57.15	Ø2.25"	~43.0	~1.69"	52.1	2.05"
CST*MW	1/16"	4, 6, 8, 10, 12, 16	Ø41.3	Ø1.625"	~35.0	~1.38"	51.4	2.025"
ST*MW	1/8"	4, 6, 8	Ø50.8	Ø2.0"	~40.0	~1.58"	52.1	2.05"
ST*MW	1/8"	10	Ø57.15	Ø2.25"	~43.5	~1.71"	52.1	2.05"
ST*MW	1/8"	12, 16	Ø61.0	Ø2.4"	~51.0	~2.00"	57.2	2.25"

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material: MATERIAL  
Dimension: DIMENSIONS

Änderungen: a) 25.02.11 FP #652, added MTG-holes (4x) ECN 15894  
MODIFICATIONS:

Projection: PROJECTION  
Allgemein: Toleranzen  
GENERAL TOLERANCES

PART NO.:

Zeichnungs-Nr. / DRAWING NO. SN EN 22768-1 / ISO 2768

Universal actuator with Multiposition valve, MW type

Masstab: 1.5:1  
Gezeichnet: 27.10.2010  
Gezeichnet: 25.02.2011  
Gezeichnet: 506384.idw

Unterwerk: KUNDE  
SUB-DIRECT:

VICI® VICI AG International, CH-6214 Schenkon  
1 - GE - 5063.84 a