



SHOWN ZRU2.5LFPK

Part#	Description	Bore		L		X1		X2		X3		X4	
		mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch	mm	Inch
ZRU1.5FPK	Union 1/16 to 1/32", fingertight	Ø0.25	Ø.010"	-45.7	-1.80"	7.9	.312"	Ø9.65	Ø.380"	Ø9.53	Ø.375"	Ø6.35	Ø.250"
ZRU1.5LFPK	Union 1/16 to 1/32", fingertight	Ø0.50	Ø.020"	-45.7	-1.80"	7.9	.312"	Ø9.65	Ø.380"	Ø9.53	Ø.375"	Ø6.35	Ø.250"
ZRU1.5TFPK	Union 1/16 to 1/32", fingertight	Ø0.79	Ø1/32"	-45.0	-1.77"	7.9	.312"	Ø9.65	Ø.380"	Ø9.53	Ø.375"	Ø6.35	Ø.250"
ZRU1.5PK	Union 1/16 to 1/32", hex	Ø0.25	Ø.010"	-35.0	-1.38"	7.9	.312"	6.35 Hex	1/4" Hex	Ø9.53	Ø.375"	Ø6.35	Ø.250"
ZRU1.5LPK	Union 1/16 to 1/32", hex	Ø0.50	Ø.020"	-35.0	-1.38"	7.9	.312"	6.35 Hex	1/4" Hex	Ø9.53	Ø.375"	Ø6.35	Ø.250"
ZRU2.5FPK	Union 1/8 to 1/32", fingertight	Ø0.25	Ø.010"	-42.6	-1.68"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	Ø6.35	Ø.250"
ZRU2.5LFPK	Union 1/8 to 1/32", fingertight	Ø0.50	Ø.020"	-42.6	-1.68"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	Ø6.35	Ø.250"
ZRU2.5TFPK	Union 1/8 to 1/32", fingertight	Ø0.79	Ø1/32"	-41.4	-1.63"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	Ø6.35	Ø.250"
ZRU21CFPK	Union 1/8 to 1/16", fingertight	Ø0.25	Ø.010"	-56.0	-2.20"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	Ø9.65	Ø.380"
ZRU21FPK	Union 1/8 to 1/16", fingertight	Ø0.75	Ø.030"	-56.0	-2.20"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	Ø9.65	Ø.380"
ZRU21LFPK	Union 1/8 to 1/16", fingertight	Ø1.00	Ø.040"	-56.0	-2.20"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	Ø9.65	Ø.380"
ZRU21TFPK	Union 1/8 to 1/16", fingertight	Ø1.59	Ø1/16"	-54.6	-2.15"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	Ø9.65	Ø.380"
ZRU21CPK	Union 1/8 to 1/16", hex	Ø0.25	Ø.010"	-45.0	-1.77"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	6.35 Hex	1/4" Hex
ZRU21MPK	Union 1/8 to 1/16", hex	Ø0.50	Ø.020"	-45.0	-1.77"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	6.35 Hex	1/4" Hex
ZRU21PK	Union 1/8 to 1/16", hex	Ø0.75	Ø.030"	-45.0	-1.77"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	6.35 Hex	1/4" Hex
ZRU21LPK	Union 1/8 to 1/16", hex	Ø1.00	Ø.040"	-45.0	-1.77"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	6.35 Hex	1/4" Hex
ZRU21TPK	Union 1/8 to 1/16", hex	Ø1.59	Ø1/16"	-43.7	-1.72"	11.1	.437"	9.53 Hex	3/8" Hex	Ø12.7	Ø.500"	6.35 Hex	1/4" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material: PEEK		Dimensional tolerances: EN ISO 2768-VS	
Revision: 5:1		Drawing No.: 04.03.2008	
Author: KUNDE		Date: 04.03.2008	
Title: Internal reducing unions high pressure PEEK		Drawing No.: 506329.idw	
Customer: KUNDE		Drawing No.: 1 - GE - 5063.29	

VICI VICI AG International, CH-6214 Schenkon