


Dimension in mm [Inch]



Prod No.:	Description	L	X
PR50A15Z1	Valco internal fittings, 1/16", 0-15 psi	~107 [-4.21"]	Hex. 6.35 [1/4"]
PR50A30Z1	Valco internal fittings, 1/16", 0-30 psi	~107 [-4.21"]	Hex. 6.35 [1/4"]
PR50A60Z1	Valco internal fittings, 1/16", 0-60 psi	~107 [-4.21"]	Hex. 6.35 [1/4"]
PR50A15Z2	Valco internal fittings, 1/8", 0-15 psi	~109.5 [-4.31"]	Hex. 9.53 [3/8"]
PR50A30Z2	Valco internal fittings, 1/8", 0-30 psi	~109.5 [-4.31"]	Hex. 9.53 [3/8"]
PR50A60Z2	Valco internal fittings, 1/8", 0-60 psi	~109.5 [-4.31"]	Hex. 9.53 [3/8"]
PR50A15E2	External fittings, 1/8", 0-15 psi	~121 [-4.76"]	Hex. 11.1 [7/16"]
PR50A30E2	External fittings, 1/8", 0-30 psi	~121 [-4.76"]	Hex. 11.1 [7/16"]
PR50A60E2	External fittings, 1/8", 0-60 psi	~123 [-4.84"]	Hex. 11.1 [7/16"]



SHOWN PR50A60Z1

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: MODIFICATIONS:		Material MATERIAL	Dimensionen DIMENSIONS
PART NO.: PR50A*		Projektion PROJECTION	Allgemein Toleranzen GENERAL TOLERANCES
			SN EN 22768-1/ISO 2768
Combo Pressure Regulators		Masstab SCALE	Gezeichnet DRAWN
		1.5:1	26.05.2009 F.P.
			Gepflicht / Freigabe CHECKED/RELEASE
			26.05.2009 FCo
			CAD Nr. CAD NO.
			506349.idw
General dimensions		Unterverz.: SUB DIREC.: KUNDE	
		Zeichnungs Nr. / DRAWING NO.	
VICI AG International, CH-6214 Schenkon		3 - GE - 5063.49	