



Ø3.6 [Ø.14"]  
8.97 [.353"]

19.05 [.75"]  
26.7 [1.05"]

Dimension in mm [Inch]



25.4 [1"]



25.15 [.99"]  
30.5 [1.2"]  
12.7 [.5"]  
19.05 [.75"]  
Ø9.53 [Ø .375"]  
Ø19.05 [Ø.75"]

Fitting size	Bore	Prod No
1/16"	0.50 mm	ASFVO
1/16"	0.75 mm	ASFVOL

Air actuated with 1" standoff

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: MODIFICATIONS:

PART NO.: ASFVO\*

Gezeichnet DRAWN: 24.02.2010 F.P.  
Geprüft / Freigabe CHECKED/ RELEASE: 24.02.2010 FCo  
CAD Nr. CAD NO.: 506358.idw

Unterverz.: SUB DIREC.: KUNDE  
Zeichnungs Nr. / DRAWING NO.: 3 - GE - 5063.58

VICI® VICI AG International, CH-6214 Schenkon