

Part#	Description	L1	L2	X1	X2
IZR1.5*	Reducer 1/16" x 1/32"	15.7 mm [.62"]	~23 mm [~.9"]	3/16" Hex	1/4" Hex
IZR21*	Reducer 1/8" x 1/16"	19.05 mm [.75"]	~28 mm [~1.1"]	1/4" Hex	3/8" Hex
IZR41*	Reducer 1/4" x 1/16"	20.8 mm [.82"]	~30 mm [~1.18"]	1/4" Hex	1/2" Hex
IZR42*	Reducer 1/4" x 1/8"	28.3 mm [1.115"]	~41 mm [~1.61"]	3/8" Hex	9/16" Hex

* different bores / different materials

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material
MATERIAL


Dimensionen
DIMENSIONS

Aenderungen:
MODIFICATIONS:

Projektion
PROJECTION

Allgemein Toleranzen
GENERAL TOLERANCES

PART NO.:



SN EN 22768-1/ ISO 2768

Internal Reducers

Masstab
SCALE

Gezeichnet
DRAWN

10.05.2005

F.P.

3:1

Geprüft / Freigabe
CHECKED/ RELEASE

12.05.2005

FCo

CAD Nr.
CAD NO.

506048.idw

General dimensions

Unterverz.:
SUB DIREC.:

KUNDE

Zeichnungs Nr. / DRAWING NO.

VICI[®]

VICI AG International, CH-6214 Schenkon

4 - GE - 5060.48