


SHOWN ZBRU21



| Part# | Description | Bore | | L | | L2 | | L3 | | L4 | | D | | X1 | X2 | X3 | X4 |
|----------|-------------------------------|-------|--------|-------|---------|-------|--------|-------|-------|------|-------|-------|---------|-----------|----------|-----------|----------|
| | | mm | Inch | mm | Inch | mm | Inch | mm | Inch | mm | Inch | mm | Inch | | | | |
| ZBRU1.5 | 1/16" to 1/32" reducing union | ∅0.25 | ∅.010" | -31.9 | -1.255" | 19.56 | .770" | 13.21 | .520" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 1/4" Hex | 3/16" Hex | 1/2" Hex |
| ZBRU1.5L | 1/16" to 1/32" reducing union | ∅0.50 | ∅.020" | -32.4 | -1.275" | 20.07 | .790" | 13.72 | .540" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 1/4" Hex | 3/16" Hex | 1/2" Hex |
| ZBRU1.5T | 1/16" to 1/32" reducing union | ∅0.89 | ∅.035" | -31.1 | -1.225" | 18.80 | .740" | 12.45 | .490" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 1/4" Hex | 3/16" Hex | 1/2" Hex |
| ZBRU2.5 | 1/8" to 1/32" reducing union | ∅0.25 | ∅.010" | -37.7 | -1.485" | 23.75 | .935" | 14.22 | .560" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 3/8" Hex | 3/16" Hex | 1/2" Hex |
| ZBRU2.5L | 1/8" to 1/32" reducing union | ∅0.50 | ∅.020" | -38.2 | -1.505" | 24.26 | .955" | 14.73 | .580" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 3/8" Hex | 3/16" Hex | 1/2" Hex |
| ZBRU2.5T | 1/8" to 1/32" reducing union | ∅0.89 | ∅.035" | -37.0 | -1.455" | 22.99 | .905" | 13.46 | .530" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 3/8" Hex | 3/16" Hex | 1/2" Hex |
| ZBRU21C | 1/8" to 1/16" reducing union | ∅0.25 | ∅.010" | -43.0 | -1.690" | 25.78 | 1.015" | 16.26 | .640" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 3/8" Hex | 1/4" Hex | 1/2" Hex |
| ZBRU21 | 1/8" to 1/16" reducing union | ∅0.75 | ∅.030" | -44.2 | -1.740" | 26.92 | 1.060" | 17.40 | .685" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 3/8" Hex | 1/4" Hex | 1/2" Hex |
| ZBRU21T | 1/8" to 1/16" reducing union | ∅1.65 | ∅.065" | -42.2 | -1.660" | 25.02 | .985" | 15.49 | .610" | 3.18 | .125" | ∅8 | ∅.3125" | 7/16" Hex | 3/8" Hex | 1/4" Hex | 1/2" Hex |
| ZBRU41C | 1/4" to 1/16" reducing union | ∅0.25 | ∅.010" | -48.4 | -1.905" | 29.59 | 1.165" | 16.89 | .665" | 4.75 | .187" | ∅11.1 | ∅.4375" | 9/16" Hex | 1/2" Hex | 1/4" Hex | 5/8" Hex |
| ZBRU41 | 1/4" to 1/16" reducing union | ∅0.75 | ∅.030" | -49.5 | -1.950" | 30.73 | 1.210" | 18.03 | .710" | 4.75 | .187" | ∅11.1 | ∅.4375" | 9/16" Hex | 1/2" Hex | 1/4" Hex | 5/8" Hex |
| ZBRU41T | 1/4" to 1/16" reducing union | ∅1.65 | ∅.065" | -47.6 | -1.875" | 28.83 | 1.135" | 16.13 | .635" | 4.75 | .187" | ∅11.1 | ∅.4375" | 9/16" Hex | 1/2" Hex | 1/4" Hex | 5/8" Hex |
| ZBRU42 | 1/4" to 1/8" reducing union | ∅0.75 | ∅.030" | -55.3 | -2.175" | 34.93 | 1.375" | 22.23 | .875" | 4.75 | .187" | ∅11.1 | ∅.4375" | 9/16" Hex | 1/2" Hex | 3/8" Hex | 5/8" Hex |
| ZBRU42L | 1/4" to 1/8" reducing union | ∅2.0 | ∅.080" | -55.3 | -2.175" | 34.93 | 1.375" | 22.23 | .875" | 4.75 | .187" | ∅11.1 | ∅.4375" | 9/16" Hex | 1/2" Hex | 3/8" Hex | 5/8" Hex |
| ZBRU42T | 1/4" to 1/8" reducing union | ∅3.25 | ∅.128" | -53.4 | -2.100" | 33.10 | 1.303" | 20.40 | .803" | 4.75 | .187" | ∅11.1 | ∅.4375" | 9/16" Hex | 1/2" Hex | 3/8" Hex | 5/8" Hex |

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen / Modifications: a) 20-01-09 FP #568, added information

| | | |
|----------|------------|---------------------|
| Material | Dimensions | EN 62278-1 ISO 2768 |
| Material | Dimensions | EN 62278-1 ISO 2768 |
| Material | Dimensions | EN 62278-1 ISO 2768 |
| Material | Dimensions | EN 62278-1 ISO 2768 |

Part No.:

Bulkhead internal reducing unions. Stainless steel

General dimensions

Material: 5.1
Scale: 5:1
Drawing No.: 20.01.2009
CAD No.: 506078.idw
User: KUNDE
Drawing No.: 1 - GE - 5060.78 a

VICI® VICI AG International, CH-6214 Schenkon