



SHOWN ECEF211.0

Part#	Description	Bore		L		X1	X2	X3
		mm	Inch	mm	Inch			
ECEF111.0*	External column end fitting 1/16" to 1/16"	Ø0.25	Ø.010"	-28.7	-1.13"	1/4" Hex	1/4" Hex	1/4" Hex
ECEF211.0*	External column end fitting 1/8" to 1/16"	Ø0.25	Ø.010"	-32.0	-1.26"	7/16" Hex	7/16" Hex	1/4" Hex
ECEF412.1*	External column end fitting 1/4" to 1/16"	Ø0.40	Ø.016"	-36.5	-1.44"	9/16" Hex	9/16" Hex	1/4" Hex
ECEF413.0*	External column end fitting 1/4" to 1/16"	Ø0.40	Ø.016"	-36.5	-1.44"	9/16" Hex	9/16" Hex	1/4" Hex
ECEF414.0*	External column end fitting 1/4" to 1/16"	Ø0.40	Ø.016"	-36.8	-1.45"	9/16" Hex	9/16" Hex	1/4" Hex
ECEF414.6*	External column end fitting 1/4" to 1/16"	Ø0.40	Ø.016"	-36.8	-1.45"	9/16" Hex	9/16" Hex	1/4" Hex
ECEF616.0*	External column end fitting 3/8" to 1/16"	Ø0.35	Ø.0145"	-36.8	-1.45"	5/8" Hex	11/16" Hex	1/4" Hex
ECEF617.0*	External column end fitting 3/8" to 1/16"	Ø0.35	Ø.0145"	-36.8	-1.45"	5/8" Hex	11/16" Hex	1/4" Hex
ECEF819.0*	External column end fitting 1/2" to 1/16"	Ø0.75	Ø.030"	-44.5	-1.75"	13/16" Hex	7/8" Hex	1/4" Hex
ECEF8110.0*	External column end fitting 1/2" to 1/16"	Ø0.75	Ø.030"	-44.5	-1.75"	13/16" Hex	7/8" Hex	1/4" Hex
ECEF1K1*	External column end fitting 1" to 1/16"	Ø0.75	Ø.030"	-53.0	-2.08"	1.25"	1 1/2" Hex	1/4" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen / MODIFICATIONS: PART NO. ECEF*		Material / MATERIAL: EN EN 22788-1 ISO 2788	Dimensional Tolerances / DIMENSIONAL TOLERANCES: F.P.
Maßstab / SCALE: 5:1	Zeichnungsnummer / Drawing No.: 21.01.2008	Stückliste / BOM No.: 506317.idw	Fertigung / PRODUCTION: F.C.o
General dimensions KUNDE		Zeichnungs Nr. / DRAWING NO.: 1 - GE - 5063.17	

VICI VICI AG International, CH-6214 Schenkon