



Dimension in mm [Inch]

Prod No	Description	Bore	L	L1
ASFV2HT	Air actuated with 2" standoff	1.50 mm	~138 [~5.44"]	50.55 [1.99"]
ASFV2HT4	Air actuated with 4" standoff	1.50 mm	~189 [~7.44"]	101.35 [3.99"]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: MODIFICATIONS: a) 24.02.10 FP #609, added ASFV2HT4
 b) 22.05.12 FP #689, changed description to 2" standoff
 c) 03.04.13 FP #717, added mounting holes ECN #18499

PART NO.: ASFV2HT*

Projektion PROJECTION:

Allgemein Toleranzen GENERAL TOLERANCES: SN EN 22768-1/ ISO 2768

Masstab SCALE: 1:1	Gezeichnet DRAWN: 10.01.2007	F.P.
	Geprüft / Freigabe CHECKED/ RELEASE: 03.04.2013	FCo
	CAD Nr. CAD NO.: 506097.idw	

General dimensions

Unterverz.: SUB DIREC.: KUNDE

Zeichnungs Nr. / DRAWING NO.: 3 - GE - 5060.97 c

VICI® VICI AG International, CH-6214 Schenkon