



Air actuated with 2" standoff <i>Prod No</i>	Air actuated with 4" standoff <i>Prod No</i>	Fitting size	Bore
ASFVOHT	ASFVOHT4	1/16"	0.50 mm
ASFVOLHT	ASFVOLHT4	1/16"	0.75 mm
ASFVO2HT	ASFVO2HT4	1/8"	1.50 mm

<i>Prod No</i>	L	L1	X
ASFVOHT	~130.5 [~5.14"]	50.55 [1.99"]	Ø31.75 [Ø1.25"]
ASFVOLHT	~130.5 [~5.14"]	50.55 [1.99"]	Ø31.75 [Ø1.25"]
ASFVO2HT	~138.0 [~5.44"]	50.55 [1.99"]	Ø36.83 [Ø1.45"]
ASFVOHT4	~181.5 [~7.14"]	101.35 [3.99"]	Ø31.75 [Ø1.25"]
ASFVOLHT4	~181.5 [~7.14"]	101.35 [3.99"]	Ø31.75 [Ø1.25"]
ASFVO2HT4	~189.0 [~7.44"]	101.35 [3.99"]	Ø36.83 [Ø1.45"]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Aenderungen: a) 22.05.12 FP #689, changed description to 2" standoff
 MODIFICATIONS: b) 15.02.13 FP #711, added ASFVO2HT*

PART NO.: ASFVO*

Gezeichnet: 25.02.2010 F.P.
 Geprüft / Freigabe: 15.02.2013 FCo
 CAD Nr.: 506362.idw

Material: DIMENSIONEN
 Allgemein Toleranzen: SN EN 22768-1/ ISO 2768

On / off valve, air actuated high temperature

Masstab SCALE: 1:1

Unterverz.: KUNDE
 SUB DIREC.: Zeichnungs Nr. / DRAWING NO.

3 - GE - 5063.62 b

VICI® VICI AG International, CH-6214 Schenkon

Dimension in mm [Inch]