


mpos. 4 port C25 valve shown

Serie	Description	Fitting	X	L
C25Z mpos.	4, 6, 8 and 10 position	1/16" ZDV	Ø34.93 mm Ø1.375"	40.0 mm 1.575"
C25Z mpos.	12 position	1/16" ZDV	Ø44.2 mm Ø1.74"	43.0 mm 1.69"
C25Z mpos.	14 position	1/16" ZDV	Ø41.3 mm Ø1.625"	40.4 mm 1.59"
C25 mpos.	4 and 6 position	1/4"-28	Ø34.93 mm Ø1.375"	40.0 mm 1.575"
C25 mpos.	8 and 10 position	1/4"-28	Ø44.2 mm Ø1.74"	40.64 mm 1.6"



THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICO. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICO AND SHALL BE RETURNED UPON DEMAND.

Material:   
 Dimensions:   
 Projection:   
 General Tolerances:   
 SN EN 22768-1 / ISO 2768

Modifications: a) 20.04.15 FP #643, New drawing

Part No.   
 Scale: 2:1   
 Drawing No. 506421.idw   
 User: KUNDE   
 Drawing No. 1 - GE - 5064.21 a

**C25 mpos valve with universal actuator**

VICO VICO AG International, CH-6214 Schenkon