


**SHOWN ZU2**



Part#	Description	Bore		L1		L2		L3		X1	X2
		mm	Inch	mm	Inch	mm	Inch	mm	Inch		
ZU.5XC	1/32" union	Ø0.15	Ø.006"	-26.0	-1.024"	17.53	.690"	-4.2	-0.167"	3/16" Hex	3/16" Hex
ZU.5	1/32" union	Ø0.25	Ø.010"	-26.0	-1.024"	17.53	.690"	-4.2	-0.167"	3/16" Hex	3/16" Hex
ZU.5L	1/32" union	Ø0.50	Ø.020"	-26.6	-1.049"	18.16	.715"	-4.2	-0.167"	3/16" Hex	3/16" Hex
ZU.5T	1/32" union	Ø0.89	Ø.035"	-25.2	-.994"	16.76	.660"	-4.2	-0.167"	3/16" Hex	3/16" Hex
ZU1XC	1/16" union	Ø0.15	Ø.006"	-37.5	-1.476"	21.59	.850"	-8.0	-0.313"	1/4" Hex	1/4" Hex
ZU1C	1/16" union	Ø0.25	Ø.010"	-37.5	-1.476"	21.59	.850"	-8.0	-0.313"	1/4" Hex	1/4" Hex
ZU1M	1/16" union	Ø0.50	Ø.020"	-37.5	-1.476"	21.59	.850"	-8.0	-0.313"	1/4" Hex	1/4" Hex
ZU1	1/16" union	Ø0.75	Ø.030"	-38.6	-1.521"	22.73	.895"	-8.0	-0.313"	1/4" Hex	1/4" Hex
ZU1L	1/16" union	Ø1.0	Ø.040"	-38.6	-1.521"	22.73	.895"	-8.0	-0.313"	1/4" Hex	1/4" Hex
ZU1T	1/16" union	Ø1.65	Ø.065"	-36.7	-1.446"	20.83	.820"	-8.0	-0.313"	1/4" Hex	1/4" Hex
ZU2C	1/8" union	Ø0.25	Ø.010"	-49.5	-1.951"	30.10	1.185"	-9.7	-0.383"	7/16" Hex	3/8" Hex
ZU2	1/8" union	Ø0.75	Ø.030"	-50.6	-1.991"	31.12	1.225"	-9.7	-0.383"	7/16" Hex	3/8" Hex
ZU2L	1/8" union	Ø2.0	Ø.080"	-50.6	-1.991"	31.12	1.225"	-9.7	-0.383"	7/16" Hex	3/8" Hex
ZU2T	1/8" union	Ø3.25	Ø.128"	-48.7	-1.916"	29.21	1.150"	-9.7	-0.383"	7/16" Hex	3/8" Hex
ZU4	1/4" union	Ø0.75	Ø.030"	-60.3	-2.375"	38.70	1.525"	-10.8	-0.425"	9/16" Hex	1/2" Hex
ZU4L	1/4" union	Ø4.6	Ø.180"	-50.6	-1.991"	31.12	1.225"	-10.8	-0.425"	9/16" Hex	1/2" Hex
ZU4T	1/4" union	Ø6.45	Ø.254"	-48.7	-1.916"	29.21	1.150"	-10.8	-0.425"	9/16" Hex	1/2" Hex
ZU6	3/8" union	Ø1.0	Ø.040"	-66.7	-2.625"	45.09	1.775"	-10.8	-0.425"	3/4" Hex	11/16" Hex
ZU6L	3/8" union	Ø7.0	Ø.276"	-66.7	-2.625"	45.09	1.775"	-10.8	-0.425"	3/4" Hex	11/16" Hex
ZU6T	3/8" union	Ø9.6	Ø.378"	-64.8	-2.550"	43.18	1.700"	-10.8	-0.425"	3/4" Hex	11/16" Hex
ZU8	1/2" union	Ø1.0	Ø.040"	-71.3	-2.805"	47.63	1.875"	-11.8	-0.465"	7/8" Hex	13/16" Hex
ZU8L	1/2" union	Ø10.0	Ø.394"	-71.3	-2.805"	47.63	1.875"	-11.8	-0.465"	7/8" Hex	13/16" Hex
ZU8T	1/2" union	Ø12.8	Ø.504"	-69.3	-2.730"	45.72	1.800"	-11.8	-0.465"	7/8" Hex	13/16" Hex

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Part No.		Material		Dimensional Tolerances	
Modifications		Production		General Tolerances	
Scale		Drawing Date		Drawing No.	
5:1		12.12.2005		F.P.	
Drawing No.		Drawing Date		Drawing No.	
506054		14.12.2005		F.C.	
Drawing No.		Drawing Date		Drawing No.	
KUNDE		506054.idw		KUNDE	
Drawing No.		Drawing Date		Drawing No.	
1 - GE - 5060.54		12.12.2005		506054.idw	

General dimensions  
VICI® VICI AG International, CH-6214 Schenkon