



SHOWN EZBRU1.5

| Part#     | Description  | Bore  |        | L     |        | L2    |       | L3    |       | L4   |       | D     |         | X1        | X2        | X3        | X4       |
|-----------|--|-------|--------|-------|--------|-------|-------|-------|-------|------|-------|-------|---------|-----------|-----------|-----------|----------|
|           |  | mm    | Inch   | mm    | Inch   | mm    | Inch  | mm    | Inch  | mm   | Inch  | mm    | Inch    |           |           |           |          |
| EZBRU1.5  | Bulkhead external / internal reducing union 1/16" to 1/32" | Ø0.25 | Ø.010" | -24.0 | - .95" | 9.78  | .385" | 5.72  | .225" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 1/4" Hex  | 3/16" Hex | 1/2" Hex |
| EZBRU1.5L | Bulkhead external / internal reducing union 1/16" to 1/32" | Ø0.50 | Ø.020" | -24.0 | - .95" | 9.78  | .385" | 5.72  | .225" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 1/4" Hex  | 3/16" Hex | 1/2" Hex |
| EZBRU1.5T | Bulkhead external / internal reducing union 1/16" to 1/32" | Ø0.89 | Ø.035" | -32.0 | -1.26" | 14.86 | .585" | 10.16 | .400" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 1/4" Hex  | 3/16" Hex | 1/2" Hex |
| EZBRU2.5L | Bulkhead external / internal reducing union 1/8" to 1/32"  | Ø0.50 | Ø.020" | -31.0 | -1.22" | 10.41 | .410" | 5.72  | .225" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 7/16" Hex | 1/4" Hex  | 1/2" Hex |
| EZBRU2.5T | Bulkhead external / internal reducing union 1/8" to 1/32"  | Ø0.89 | Ø.035" | -31.0 | -1.22" | 10.41 | .410" | 5.72  | .225" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 7/16" Hex | 1/4" Hex  | 1/2" Hex |
| EZBRU21   | Bulkhead external / internal reducing union 1/8" to 1/16"  | Ø0.75 | Ø.030" | -37.5 | -1.48" | 13.72 | .540" | 9.02  | .355" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 7/16" Hex | 1/4" Hex  | 1/2" Hex |
| EZBRU21L  | Bulkhead external / internal reducing union 1/8" to 1/16"  | Ø1.0  | Ø.040" | -37.5 | -1.48" | 13.72 | .540" | 9.02  | .355" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 7/16" Hex | 1/4" Hex  | 1/2" Hex |
| EZBRU21T  | Bulkhead external / internal reducing union 1/8" to 1/16"  | Ø1.65 | Ø.065" | -37.5 | -1.48" | 13.72 | .540" | 9.02  | .355" | 3.18 | .125" | Ø8.0  | Ø.3125" | 7/16" Hex | 7/16" Hex | 1/4" Hex  | 1/2" Hex |
| EZBRU41   | Bulkhead external / internal reducing union 1/4" to 1/16"  | Ø0.75 | Ø.030" | -45.5 | -1.79" | 18.03 | .710" | 12.37 | .487" | 4.75 | .187" | Ø11.1 | Ø.4375" | 9/16" Hex | 9/16" Hex | 1/4" Hex  | 5/8" Hex |
| EZBRU41T  | Bulkhead external / internal reducing union 1/4" to 1/16"  | Ø1.65 | Ø.065" | -45.5 | -1.79" | 18.03 | .710" | 12.37 | .487" | 4.75 | .187" | Ø11.1 | Ø.4375" | 9/16" Hex | 9/16" Hex | 1/4" Hex  | 5/8" Hex |
| EZBRU42   | Bulkhead external / internal reducing union 1/4" to 1/8"   | Ø0.75 | Ø.030" | -47.2 | -1.86" | 18.03 | .710" | 12.37 | .487" | 4.75 | .187" | Ø11.1 | Ø.4375" | 9/16" Hex | 9/16" Hex | 3/8" Hex  | 5/8" Hex |
| EZBRU42L  | Bulkhead external / internal reducing union 1/4" to 1/8"   | Ø2.0  | Ø.080" | -47.2 | -1.86" | 18.03 | .710" | 12.37 | .487" | 4.75 | .187" | Ø11.1 | Ø.4375" | 9/16" Hex | 9/16" Hex | 3/8" Hex  | 5/8" Hex |
| EZBRU42T  | Bulkhead external / internal reducing union 1/4" to 1/8"   | Ø3.18 | Ø.125" | -47.2 | -1.86" | 18.03 | .710" | 12.37 | .487" | 4.75 | .187" | Ø11.1 | Ø.4375" | 9/16" Hex | 9/16" Hex | 3/8" Hex  | 5/8" Hex |

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REPRODUCED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen / Modifications: a) 19.01.09 FP #568, added information

Prohibition of Production: KUNDE

Part No.:

Material:

Dimensions:

Scale: 7:1

Customer Order No.: 10.03.2006 F.P.

Production Date: 19.01.2009 F.C.o

Customer Order No.: 506075.idw

Customer: KUNDE

Drawing No. / Drawing No.: 1 - GE - 5060.75 a

VICI VICI AG International, CH-6214 Schenkon