


Shown EAOR21



Shown EAOR22

Part#	Description	Bore		L		X1	X2
		mm	Inch	mm	Inch		
EAOR21	External adaptor 1/8" to 10-32 thread	Ø1.58	Ø.062"	~25.9	~1.02	7/16" Hex	3/8" Hex
EAOR22	External adaptor 1/8" to 5/16-24 thread	Ø1.58	Ø.062"	~27.4	~1.08	7/16" Hex	
EAOR22L	External adaptor 1/8" to 5/16-24 thread	Ø2.16	Ø.085"	~27.4	~1.08	7/16" Hex	

Dimension in mm [Inch]

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Änderungen:
MODIFICATIONS:
PART NO.:

Material	Dimensions
PROJECTION	GENERAL TOLERANCES
SN EN 22768-1/ISO 2768	

**Adaptor
1/8" external to ***

General dimensions

Maßstab SCALE	Gezeichnet DRAWN	22.06.2006	F.P.
4:1	Geprüft / Freigegeben CHECKED / RELEASED	22.06.2006	FCo
	CAD Nr. CAD NO.	506084.idw	

Unterwerk:
SUB DIR.: **KUNDE**

VICI® VICI AG International, CH-6214 Schenkon

Zeichnungs-Nr. / DRAWING NO.
2 - GE - 5060.84