






	Valve Type	Connections	No. of Pos.	X		Y		L	
				mm	Inch	mm	Inch	mm	Inch
 <p>ONE OFFSET PORT</p> <p>ONE ROW OF PORTS</p>	CSD*MW	1/16"	4, 6, 8, 10, 12, 16	ø38.1	ø1.5	~34.5	~1.36	41.3	1.625
	SD*MW	1/8"	4, 6, 8, 10, 12	ø50.8	ø2	~40.5	~1.60	41.9	1.65
	SD*MW	1/8"	16	ø57.15	ø2.25	~51.0	~2.00	41.9	1.65
 <p>TWO OFFSET PORTS</p> <p>ONE ROW OF PORTS</p>	CSC*MW	1/16"	4, 6, 8, 10, 12, 16	ø41.3	ø1.625	~36.0	~1.42	51.4	2.025
	SC*MW	1/8"	4, 6, 8, 10	ø50.8	ø2	~38.5	~1.52	52.1	2.05
	SC*MW	1/8"	12, 16	ø57.15	ø2.25	~43.0	~1.69	52.1	2.05
 <p>ONE OFFSET PORT</p> <p>TWO ROWS OF PORTS</p>	CSF*MW	1/16"	4, 6, 8, 10, 12, 16	ø41.3	ø1.625	~35.0	~1.38	51.4	2.025
	SF*MW	1/8"	4, 6, 8, 10	ø50.8	ø2	~39.5	~1.56	52.1	2.05
	SF*MW	1/8"	12, 16	ø57.15	ø2.25	~43.0	~1.69	52.1	2.05
 <p>TWO OFFSET PORTS</p> <p>TWO ROWS OF PORTS</p>	CST*MW	1/16"	4, 6, 8, 10, 12, 16	ø41.3	ø1.625	~35.0	~1.38	51.4	2.025
	ST*MW	1/8"	4, 6, 8	ø50.8	ø2	~40.0	~1.58	52.1	2.05
	ST*MW	1/8"	10	ø57.15	ø2.25	~43.5	~1.71	52.1	2.05
	ST*MW	1/8"	12, 16	ø61	ø2.4	~51.0	~2.00	57.2	2.25



Dimension in mm [Inch]

Änderungen: a) 07.09.09 FP #592, chg. table data Y	
PART NO.:	
<b>Micro electric actuator</b>	19.05.1999 F.P.
<b>Multiposition Valve</b>	1.5:1 07.09.2009 FCo
„Dimensions	508024.DWG
Sub title: KUNDE	
Zeichnung Nr. / DRAWING NO.	
VICI	1 GE 5060 24 a

VICI AG International, CH-6214 Schenkon