



Dimension in mm [Inch]



SHOWN EUTA-CSC6MWE

Type	Connections	No. of Pos.	X		Y		L	
			mm	Inch	mm	Inch	mm	Inch
CSD*MW	1/16"	4, 6, 8, 10, 12, 16	Ø38.1	Ø1.5"	~34.5	~1.36"	41.3	1.625"
SD*MW	1/8"	4, 6, 8, 10, 12	Ø50.8	Ø2.0"	~40.5	~1.60"	41.9	1.65"
SD*MW	1/8"	16	Ø57.15	Ø2.25"	~51.0	~2.00"	41.9	1.65"
CSC*MW	1/16"	4, 6, 8, 10, 12, 16	Ø41.3	Ø1.625"	~36.0	~1.42"	51.4	2.025"
SC*MW	1/8"	4, 6, 8, 10	Ø50.8	Ø2.0"	~38.5	~1.52"	52.1	2.05"
SC*MW	1/8"	12, 16	Ø57.15	Ø2.25"	~43.0	~1.69"	52.1	2.05"
CSF*MW	1/16"	4, 6, 8, 10, 12, 16	Ø41.3	Ø1.625"	~35.0	~1.38"	51.4	2.025"
SF*MW	1/8"	4, 6, 8, 10	Ø50.8	Ø2.0"	~39.5	~1.56"	52.1	2.05"
SF*MW	1/8"	12, 16	Ø57.15	Ø2.25"	~43.0	~1.69"	52.1	2.05"
CST*MW	1/16"	4, 6, 8, 10, 12, 16	Ø41.3	Ø1.625"	~35.0	~1.38"	51.4	2.025"
ST*MW	1/8"	4, 6, 8	Ø50.8	Ø2.0"	~40.0	~1.58"	52.1	2.05"
ST*MW	1/8"	10	Ø57.15	Ø2.25"	~43.5	~1.71"	52.1	2.05"
ST*MW	1/8"	12, 16	Ø61.0	Ø2.4"	~51.0	~2.00"	57.2	2.25"

THIS DOCUMENT AND THE INFORMATION WHICH IT CONTAINS SHALL NOT BE USED, EXPLOITED OR SOLD, AND SHALL NOT BE REVEALED OR DISCLOSED TO OTHERS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VICI. THIS DOCUMENT SHALL REMAIN THE PROPERTY OF VICI AND SHALL BE RETURNED UPON DEMAND.

Material: MATERIAL
Dimension: DIMENSIONS

Änderungen: a) 25.02.11 FP #652, added MTG-holes (4x) ECN 15894
MODIFICATIONS:

Projection: PROJECTION
Allgemein: Toleranzen
GENERAL TOLERANCES

PART NO.:

Scale: 1.5:1

Universal actuator with Multiposition valve, MW type

General dimensions

Customer: KUNDE

Zeichnungs Nr. / DRAWING NO. 1 - GE - 5063.84 a

VICI® VICI AG International, CH-6214 Schenkon

27.10.2010 F.P.
25.02.2011 FCo
506384.idw

SN EN 22768-1 / ISO 2768